

AQUACLEAN mini

MAKE YOUR SKIN YOUNG

Pore
Cleansing

Facial
Lifting

Solution
Delivery

With valve
finely adjustable
amount

3 The combination of
functions
for complete skin care



SPECIFICATION

Applicators	Aqua cleansing / Electroporation Double Lifting
Input Voltage	AC 100~120V/220~240V, 50/60Hz
Weight(kg)	17kg
Dimension(WxDxH mm)	430 x 380 x 380



EunSung Global Corp.

#1102~1105, 53, Gasan digital 2-ro,
Geumcheon-gu, Seoul, 153-706, Korea
Tel: +82-2-514-9713(Rep) Fax: +82-2-514-9716
e-mail: es@esglobal.co.kr/www.esglobal.co.kr

 EunSung

AQUACLEAN mini

MAKE YOUR SKIN YOUNG

3-in-1

AquaClean mini consist of 3-in-1 system that combines three different technologies including aqua clean system using vacuum, electroporation and double lifting. 3-in-1 system provides all facial care including cleansing, solution penetration and lifting to a single device.



① Aqua Cleansing



② Multi-polarity Electroporation



③ Double Lifting

Valve type
Vacuum & Flow rate control

Cleansing & Exfoliation
by vacuum



Removal effect of sebum,
white head, black head



Effective skin cleansing



Various solutions
can be Applied



The effects of specialized
facial lifting treatment



Penetration of sloution
to the deep dermis



AQUACLEAN mini 3-in-1 System

Optimum Multi-functional skin treatment

Hypoallergenic
peeling function


+

Solution
Penetration

+

Double lifting

1 Aqua Clean system



- Hypoallergenic peeling function that removal of sebum, dead skin cell and wastes using vacuum.
- Simple treatment procedure without any side-effects.
- Smooth peeling & treatment without stimulate.



Small



Large



Cap

2 Electroporation

- Without damage to the cellular skin structure, micro-current creates a temporary opening on the skin for penetration of nourishing ampoules into deep layers of the skin.
- Adjustable Treatment Level (Level 1 to Level 10)

3 Double Lifting

- By supplying micro-current to flow through superficial layer of epidermis, lifting effect from stimulated facial muscle is achieved.
- Adjustable Treatment Level (Level 1 to Level 10)


Premium Solution



1 Remove sebum and oxidized waste deeper in the pores by addition of Lactic Acid.



2 Dissolve oxidized dead skin and calms ACNE, rash, and other skin problem by addition of Salicylic Acid.



3 Supplement nutrients and deep moisturizing on the cleaner skin by addition of Portulaca leracea extract.



4 Containing ethanol, solution 4 disinfects and cleanses the entire system.